

Advanced CROW Order Form

 Left Foot **Right Foot** Note: Please complete our regular order form if you require a custom molded shoe for the other foot

 42 Breckenridge Street • Buffalo, NY 14213 • (716) 881-4484 • Fax (716) 881-0406 • info@trumold.com
COMPANY: DATE _____ PO. # _____

Name _____

Phone _____ Fax _____

E Mail _____

SHIP TO: 2-Day Ground Overnight

Name _____

Address _____

City _____ State _____ Zip _____

BILL TO:

Name _____

Address _____

City _____ State _____ Zip _____

PATIENT INFO:

PT Name: _____

Address _____

City _____ State _____ Zip _____

 Weight _____ Age _____ Sex: M F

Occupation _____

Had Tru-Mold Shoes before

(Approx. Date) _____ Invoice/CAD/CAM number _____

INDICATE ANY LESIONS:

(i.e. Amputations or Ulcerations)

Left
Right

Dorsum
Plantar
Plantar
Dorsum
DIAGNOSIS:
 Diabetic Charcot Club Foot

 Amputee (where): _____

 Other _____

STANDARD FEATURES:

- Black leather
- Regular lace eyelets from bottom to ankle then speed-lace up to topline
- Surgical opening
- Plastazote lining throughout
- One custom molded triple density insert
 - 1/4" Pink Plastazote top
 - 1/8" Poron middle
 - 1/4" White Plastazote bottom
- 1" Heel-to-toe rocker sole with non-slip final outsole
- V Factor Plate – Tru-Mold's exclusive varus/valgus auxiliary control device
(see back of order form for further details)
- TM Pressure Relief System – Tru-Mold's exclusive method to reduce pressure on ulcerations, bone protrusions and other problem areas – includes a removable poron plug for additional relief
(see back of order form for further details)

FINISHED INSIDE HEIGHT:

(Includes height of plastic from bottom of barefoot heel and 1" rolled padded collar) _____

CAST MODIFICATION:
 Reg. Elongation (5/8") 1/8" Extra (5/8" + 1/8")

 High Toe Box Extra High Toe Box

 Other _____

ANKLE RANGE OF MOTION:

 Dorsiflexion - Fixed/Fused Mobile/Normal

 Plantarflexion - Fixed/Fused Mobile/Normal

 Hindfoot - Flexible - Correct to Neutral

 Rigid - Do Not Correct / Requires Extrinsic Posting

 Forefoot - Flexible - Correct to Neutral

 Rigid - Do Not Correct / Requires Extrinsic Posting

PLACE A CHECK MARK IF YOU WISH THE FOLLOWING FREE OF CHARGE OPTIONS:

- Velcro D-ring (external) with 1” straps instead of lace/speed-lace
- Leather lining instead of Plastazote lining
- Forefoot only rocker sole instead of heel-to-toe rocker

The following options are available for an additional charge:

- | | LF | RF |
|--|-----|-----|
| <input type="checkbox"/> One extra custom molded insert | ___ | ___ |
| <input type="checkbox"/> Two extra custom molded inserts | ___ | ___ |
| <input type="checkbox"/> Medial flare | ___ | ___ |
| <input type="checkbox"/> Lateral flare | ___ | ___ |
| <input type="checkbox"/> Amputation Filler | ___ | ___ |

Short Limb Raise:

- | | LF | RF |
|-------------------------------|-----|-----|
| <input type="checkbox"/> Heel | ___ | ___ |
| <input type="checkbox"/> Ball | ___ | ___ |
| <input type="checkbox"/> Toe | ___ | ___ |

Note – All raises placed on outside with soling

Exclusive features of Tru-Mold’s Advanced CROW:

V Factor Plate – this separate device is not lasted into the Advanced CROW at the time of manufacture. It is intended to provide auxiliary varus/valgus correction if your patient’s condition changes at a later date.

To use the V Factor Plate, just lift out the insert and the removable base material, secure the V Factor Plate with double-sided tape to the desired location on the bottom of the boot and place the base material and insert back into the Advanced CROW.

TM Pressure Relief System – the Advanced CROW will be recessed on the insert and removable base material to reduce the plantar pressure on ulcerations, bone protrusions and other problem areas.

If your patient requires additional pressure relief at these sites, simply remove the poron plug from the base material.

SPECIAL “FOOT NOTES”

WEIGHT-BEARING TRACINGS

(Hold pencil vertically when tracing)

Left Heel

Right Heel